## **FEROS Thorium Argon Calibration Lamps Log**

Thorium Argon calibration lamps used in FEROS must be of high purity type (HP) type 4060AHP, small window (25 mm)

ThAr lamps list with respective S/N for each lamp received on January 2011. The spares will be kept in the store room, ground floor NOB building.

Serial	Type	Status	When	Who	Comments
No.					
		In use			
531290	4160AHP	FEROS		AGO	Used for test
531402	4160AHP	Stock			
531485	4160AHP	Stock			
531488	4160AHP	Stock			
531489	4160AHP	Stock			
531490	4160AHP	Stock			
531491	4160AHP	Stock			
531493	4160AHP	Stock			
531497	4160AHP	Stock			
531498	4160AHP	Stock			